

## Product datasheet for AM50111PU-S

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Fibronectin (FN1) Mouse Monoclonal Antibody [Clone ID: SPM539]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: SPM539
Applications: FC, IF, IHC

Recommended Dilution: Flow Cytometry:  $0.5-1 \mu g/10^6$  cells.

Immunofluorescence: 0.5-1 µg/ml.

Immunohistochemistry on Formalin-Fixed Paraffin Sections: 0.5-1 µg/ml for 30 minutes

at RT.

Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH

6.0, for 10-20 min followed by cooling at RT for 20 minutes.

**Positive Control**: SW156 cells or Kidney.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human fibronectin purified from serum by affinity chromatography on gelatin-sepharose.

**Specificity:** This MAb is directed against the peptide core and reacts with both the plasma and cellular

forms of fibronectin. It blocks the fibronectin-medicated cell attachment not by disrupting the collagen-fibronectin interaction, but by interfering with the attachment of fibronectin to its

receptor on the cell surface.

*Cellular Localization*: Connective tissue matrix.

Formulation: 10mM PBS

State: Purified

State: Liquid purified IgG fraction from Bioreactor Concentrate

Stabilizer: 0.05% BSA

Preservative: 0.05% Sodium Azide

**Concentration:** lot specific

**Purification:** Protein A/G Chromatography

Conjugation: Unconjugated

Storage: Store undiluted at 2-8°C.





**Stability:** Shelf life: one year from despatch.

Predicted Protein Size: 220 kDa (monomer); 440 kDa (dimer)

**Gene Name:** fibronectin 1

Database Link: Entrez Gene 2335 Human

P02751

**Background:** Fibronectin is an extracellular matrix glycoprotein present on most cell surfaces, in

extracellular fluids and in plasma. A high molecular weight heterodimeric protein, it was originally discovered as a protein missing from the surfaces of virus-transformed cells, and it has been shown to be involved in various functions including cell adhesion, cell motility and wound healing. Alternative splicing and glycosylation give rise to several different forms of

Fibronectin, some of which exhibit restricted tissue distribution or association with

malignancies. It has been shown that Myofibroblast phenotype formation correlates with the occurrence of glycosylated Fibronectin and Fibronectin splice variants in Dupuytren's disease. Fibronectin is a dimeric glycoprotein of 440kDa, which is present in cells, extracellular matrix,


and blood. It possesses at least four binding sites for collagen, glycosaminoglycans,

transglutaminase, and a cell surface receptor. Fibronectin is involved in cell adhesion, tissue

organization, and wound healing.

Synonyms: FN1, Cold-insoluble globulin, CIG

## **Product images:**



Formalin-fixed, paraffin-embedded Human Pancreatic Adenocarcinoma stained with Fibronectin Monoclonal Antibody (Clone SPM539).