

## Product datasheet for **AM50114PU-T**

### HSP27 (HSPB1) Mouse Monoclonal Antibody [Clone ID: SPM252]


#### Product data:

Product Type:	Primary Antibodies
Clone Name:	SPM252
Applications:	IF, IHC, IP, WB
Recommended Dilution:	<b>ELISA:</b> Use Antibody without BSA for Coating. <b>Flow Cytometry:</b> 0.5-1 µg/10 <sup>6</sup> cells. <b>Immunofluorescence:</b> 0.5-1 µg/ml. <b>Western Blot:</b> 0.5-1 µg/ml. <b>Immunoprecipitation:</b> 0.5-1 µg/500 µg protein lysate. <b>Immunohistochemistry on Frozen and Formalin-Fixed Paraffin Sections:</b> 0.5-1 µg/ml for 30 minutes at RT. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes. <b>Positive Control:</b> ~50% of Breast Carcinomas are positive for hsp27 especially those that are also positive for estrogen and/or progesterone receptor.
Reactivity:	Chimpanzee, Chicken, Human, Mouse, Rat, Sheep
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Partially purified hsp27 (earlier called 24K) protein from breast cancer MCF-7 cells.
Specificity:	It recognizes a 24-27kDa estrogen-regulated protein, identified as heat shock protein 27 (hsp27). Hsp27 was recently found to be identical to the estrogen-induced "p29" and "24K" protein. About 50% of breast carcinomas are positive for hsp27 especially those that are also positive for estrogen and/or progesterone receptor. HSP27 has also been implicated in drug resistance in cancer cells. <b>Cellular Localization:</b> Predominantly cytoplasmic with some nuclear.
Formulation:	10mM PBS State: Purified State: Liquid purified IgG fraction from Bioreactor Concentrate Stabilizer: 0.05% BSA Preservative: 0.05% Sodium Azide




[View online »](#)

Concentration:	lot specific
Purification:	Protein A/G Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	27 kDa
Gene Name:	heat shock protein family B (small) member 1
Database Link:	<a href="#">Entrez Gene 3315 Human P04792</a>
Synonyms:	Heat shock protein beta-1, Heat shock 27 kDa protein, HSP28, 28 kDa heat shock protein, SRP27, HSP25

**Product images:**

Formalin-Paraffin breast carcinoma (10X) stained with HSP27 Antibody (Clone SPM252).



Formalin-Paraffin breast carcinoma (20X) stained with HSP27 Antibody (Clone SPM252).