

Product datasheet for **AM50119PU-S**

Human IgM (heavy chain) Mouse Monoclonal Antibody [Clone ID: SPM557]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM (heavy chain) Mouse Monoclonal Antibody [Clone ID: SPM557]
Clone Name:	SPM557
Applications:	FC, IF, IHC, IP, WB
Recommended Dilution:	ELISA: Use Antibody without BSA for Coating. Flow Cytometry: 0.5-1 µg/10 ⁶ cells. Immunofluorescence: 0.5-1 µg/ml. Western Blot: 0.5-1 µg/ml. Immunoprecipitation: 0.5-1 µg/500 µg protein lysate. Immunohistochemistry on Frozen and Formalin-Fixed Paraffin Sections: 0.5-1 µg/ml for 30 minutes at RT. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 minutes. Positive Control: 293T, Raji or hPBL cells. Tonsil or Spleen.
Reactivity:	Human
Host:	Mouse
Immunogen:	Recombinant heavy chain of human IgM.
Isotype:	IgG1
Formulation:	10mM PBS State: Purified State: Liquid purified IgG fraction from Bioreactor Concentrate Stabilizer: 0.05% BSA Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Purification:	Protein A/G Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.



[View online »](#)