

Product datasheet for **AM50286PU-S**

IGHG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: ICO-97]

Product data:

Product Type:	Secondary Antibodies
Product Name:	IGHG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: ICO-97]
Clone Name:	ICO-97
Applications:	FC, IF, IHC
Recommended Dilution:	ELISA: Use Antibody without BSA for Coating. Flow Cytometry: 0.5-1 µg/10 ⁶ cells. Immunohistochemistry on Frozen Sections: 0.5-1 µg/ml for 30 minutes at RT. Positive Control: 293T, Raji or hPBL cells. Tonsil or Spleen.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human Ig Gamma Chain.
Isotype:	IgG1
Formulation:	10mM PBS State: Purified State: Liquid purified IgG fraction from Bioreactor Concentrate Stabilizer: 0.05% BSA Preservative: 0.05% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.



[View online »](#)