

## Product datasheet for AM50286PU-T

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## IGHG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: ICO-97]

**Product data:** 

**Product Type:** Secondary Antibodies

Product Name: IGHG1 (heavy chain) Mouse Monoclonal Antibody [Clone ID: ICO-97]

Clone Name: ICO-97
Applications: FC, IF, IHC

**Recommended Dilution: ELISA**: Use Antibody without BSA for Coating.

Flow Cytometry:  $0.5-1 \mu g/106$  cells.

**Immunohistochemistry on Frozen Sections:** 0.5-1 µg/ml for 30 minutes at RT.

Positive Control: 293T, Raji or hPBL cells. Tonsil or Spleen.

Reactivity: Human Host: Mouse

**Immunogen:** Purified human Ig Gamma Chain.

**Isotype:** IgG1

**Formulation:** 10mM PBS

State: Purified

State: Liquid purified IgG fraction from Bioreactor Concentrate

Stabilizer: 0.05% BSA

Preservative: 0.05% Sodium Azide

Concentration: lot specific

**Conjugation:** Unconjugated

**Storage:** Store undiluted at 2-8°C.