

## **Product datasheet for AM50624PU-N**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Carbonyl reductase 1 Mouse Monoclonal Antibody [Clone ID: AT2D6]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: AT2D6

**Applications:** ELISA, FC, IF, WB

Recommended Dilution: The antibody has been tested by ELISA, Western blot analysis, ICC/IF and Flow cytometry to

assure specificity and reactivity. Since application varies, however, each investigation should

be titrated by the reagent to obtain optimal results.

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

**Immunogen:** Recombinant human CBR1 (1-277aa) purified from E. coli.

**Specificity:** Recognizes Human CBR1.

**Formulation:** Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

State: Purified

State: Liquid purified Ig fraction

**Concentration:** lot specific

**Purification:** Protein-A affinity chromatography

**Conjugation:** Unconjugated

Storage: Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

**Gene Name:** carbonyl reductase 1

Database Link: Entrez Gene 873 Human

P16152





**Background:** Carbonyl reductase 1 (CBR1) is a DPH-dependent, monomeric, and cytosolic enzyme


belonging to a family of short-chain dehydrogeses/reductases. This protein consists of 277 amino acid residues and is widely distributed in human tissues such as liver, epidermis, stomach, small intestine, kidney, neurol cells, and smooth muscle fiber. CBR1 metabolizes many toxic environmental quinones and pharmacological relevant substrates such as the anticancer drug, doxorubicin. The best substrates of CBR1 are quinones, including ubiquinone-1 and tocophrolquinone (vitamin E).

Synonyms: CBR1, CBR, CRN


**Protein Families:** Druggable Genome

**Protein Pathways:** Arachidonic acid metabolism, Metabolic pathways

## **Product images:**



The Recombinant Human CBR1 (50ng) and Cell lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human CBR1 antibody (1:1000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: Recombinant Human CBR1 Lane 2.: HeLa cell lysate Lane 3.: 293T cell lysate Lane 4.: MCF-7 cell lysate Lane 5.: HepG2 cell lysate



ICC/IF analysis of CBR1 in HeLa cells line, stained with DAPI (Blue) for nucleus staining and monoclonal anti-human CBR1 antibody (1:100) with goat anti-mouse IgG-Alexa fluor 488 conjugate (Green).