

Product datasheet for **AM50632PU-N**

PMM2 Mouse Monoclonal Antibody [Clone ID: AT3B4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	AT3B4
Applications:	ELISA, FC, WB
Recommended Dilution:	The antibody has been tested by ELISA, Western blot analysis and Flow cytometry to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Recombinant human PMM2 (1-246aa) purified from E.coli.
Formulation:	Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol. State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein-A affinity chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	phosphomannomutase 2
Database Link:	Entrez Gene 5373 Human O15305



[View online »](#)


Background: PMM2, also known as Phosphomannomutase 2, belongs to the eukaryotic PMM family. PMM2 is involved in the synthesis of the GDP-mannose and dolichol-phosphate-mannose required for a number of critical mannosyl transfer reactions. It catalyzes the isomerization of mannose 6-phosphate to mannose 1-phosphate. Mutations in PMM2 are associated with Congenital disorders of glycosylation (CDG)-Ia, an autosomal recessive disorder characterized by central nervous system dysfunction and multiorgan failure.

Synonyms: Phosphomannomutase 2, PMM 2

Protein Families: Druggable Genome

Protein Pathways: Amino sugar and nucleotide sugar metabolism, Fructose and mannose metabolism, Metabolic pathways

Product images:



The Cell lysates (40ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human PMM2 antibody (1:500). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: A549 cell lysate
Lane 2.: HepG2 cell lysate