

Product datasheet for **AP00683PU-N**

Clostridium (genus specific) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF
Recommended Dilution:	Suitable for use in IFA and conjugation purposes.
Reactivity:	Clostridium perfringens, Clostridium sordellii, Clostridium novyi, Clostridium septicum, Clostridium chauvoei
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Extract of Cl. perfringens, Cl. sordellii, Cl. novyi, Cl. septicum and Cl. chauvoei
Specificity:	This antibody is reactive with proteins from common Clostridia species. Antiserum is not absorbed and may react with related microorganisms. Reactive (IFA) with Cl. perfringens, Cl. sordellii, Cl. novyi, Cl. septicum, Cl. chauvoei.
Formulation:	0.01M PBS, pH 7.2 containing 0,09% Sodium Azide as preservative State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Note:	Centrifuge before opening to ensure complete recovery of vial contents.



[View online »](#)