

Product datasheet for **AP00869PU-N**

Staphylococcus aureus Enterotoxin D Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Staphylococcus aureus Enterotoxin D (SED)
Specificity:	This antibody is specific to Staphylococcus aureus Enterotoxin D (SED). Minimal cross-reactivity with staphylococcal enterotoxins A, B, C, E, ET and TSST and alpha hemolysin.
Formulation:	0.01M phosphate, pH 7.4 containing 0.9% Sodium chloride without preservatives State: Aff - Purified State: Lyophilized purified Ig fraction
Reconstitution Method:	Restore with 1ml deionized water
Concentration:	lot specific
Purification:	Toxin Specific Immunoaffinity Chromatography
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antibody at -20°C. If aliquoted for long term storage, fill volume should be equal to or greater than 50% of the nominal fill volume of the vial used. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P20723
Background:	Staphylococcus aureus are gram-positive bacteria that, although part of the normal flora of humans (nasal passages, skin and mucous membranes), have the potential to become opportunistic pathogens. S. aureus is a frequent cause of a wide range of infections and some bacterial strains have become resistant to many commonly used antibiotics.
Synonyms:	SED, entD, S. aureus



[View online »](#)