

Product datasheet for **AP01109BT-S**

IGFBP7 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA: Direct: To detect hIGF-BP7 by direct ELISA (using 100 µl/well antibody solution) a concentration of 0.25 - 1.0 µg/ml is required. In conjunction with compatible secondary reagents, it allows the detection of at least 0.2 - 0.4 ng/well of recombinant hIGF-BP7. Western Blot: To detect hIGF-BP7 by Western Blot analysis this antibody can be used at a concentration of 0.1 - 0.2 µg/ml. Used in conjunction with compatible secondary reagents the detection limit for recombinant hIGF-BP7 is 1.5 - 3.0 ng/lane, under either reducing or non-reducing conditions.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Highly pure (> 98 %) recombinant human IGFBP7
Specificity:	This antibody detects IGFBP7
Formulation:	PBS, pH 7.2 Label: Biotin State: Sterile filtered, lyophilized Ig fraction
Reconstitution Method:	Centrifuge vial prior to opening. Restore in sterile PBS containing 0.1 % BSA to a concentration of 0.1 - 1.0 mg/ml.
Purification:	Affinity purified
Conjugation:	Biotin
Storage:	Store the lyophilized antibody at -20 °C. Following reconstitution it is stable for two weeks at 2 - 8 °C. Frozen aliquots are stable for 6 months when stored at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	insulin like growth factor binding protein 7
Database Link:	Entrez Gene 3490 Human Q16270



[View online »](#)

Background:

IGF-BPs form high affinity complexes with both IGF-I and IGF-II and act to control of the distribution, function and activity of IGFs in various cell tissues and body fluids. There are seven named IGF-BPs. IGFBP7 plays a role in skeletal myogenesis by binding to IGF in a manner that inhibits IGF induced differentiation of skeletal myoblasts, without affecting IGF induced proliferation. Additionally, IGFBP7 suppresses growth and colony formation of prostate and breast cancer cell lines through an IGF independent mechanism, which causes a delay in the G1 phase of the cell cycle, and increased apoptosis. IGFBP7 is expressed in a wide range of normal human tissues and it usually shows reduced expression in cancer cell lines of prostate, breast, colon, and lung origin.

Synonyms:

IGF-binding protein 7, IGFBP-7, IBP7, IBP-7, IGFBP-rP1, MAC25, PSF, TAF