

Product datasheet for AP01225PU-M

KCNJ11 Rabbit Polyclonal Antibody

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Recommended Dilution:	ELISA: 1/10000-1/20000. Western Blot: 1/500-1/1000. Immunofluorescence: 1/50-1/200. Immunohistochemistry on Paraffin Sections: 1/50-1/200.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide, corresponding to amino acids 188-243 of Human Kir6.2.
Specificity:	This antibody detects endogenous levels of Kir6.2 protein. (region surrounding Val220)
Formulation:	Phosphate buffered saline (PBS), pH~7.2 containing 0.05% Sodium Azide as preservative. State: Aff - Purified State: Liquid purified Ig fraction (> 95% pure by SDS-PAGE)
Concentration:	1.0 mg/ml
Purification:	Affinity Chromatography using epitope-specific immunogen.
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Predicted Protein Size:	~40 kDa
Gene Name:	potassium voltage-gated channel subfamily J member 11
Database Link:	<u>Entrez Gene 3767 Human</u> <u>Q14654</u>



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2025 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

GRIGENE KCNJ11 Rabbit Polyclonal Antibody – AP01225PU-M

Background:ATP-sensitive K+ channels play important roles in many cellular functions by coupling cell
metabolism to electrical activity. KIR6.1 and KIR6.2 are members of the KIR (inwardly
rectifying potassium channel) family of potassium channels. Inward rectifying K+ channels
possess a greater tendency to allow potasium to flow into the cell rather than out of it. These
channels comprise two subunits: a KIR6.0 subfamily component and a SUR component, which
is a member of the ATP-binding cassette protein superfamily. Mutations in the gene coding
for these channels are a cause of an autosomal recessive disorder characterized by
unregulated insulin secretion. The amino-terminal and carboxyl-terminal domains of KIR
channel subunits are both intracellular, and the two intracellular domains of KIR6.2 physically
interact with each other.

Synonyms:

KCNJ11, Inward rectifier K+ channel Kir6.2, IKATP

Product images:



Western blot (WB) analysis of Kir6.2 antibody (Cat.-No.: [AP01225PU-N]) in extracts from HeLa cells.



Immunohistochemistry (IHC) analyzes of Kir6.2 antibody (Cat.-No.: [AP01225PU-N]) in paraffinembedded human brain tissue.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2025 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US