

## Product datasheet for **AP01366PU-M**

### **TIS11B (ZFP36L1) Rabbit Polyclonal Antibody**

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Recommended Dilution:	<b>Western Blot:</b> 1/500-1/1000. <b>Immunohistochemistry on paraffin sections:</b> 1/50-1/200.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	TIS11B antibody detects endogenous levels of TIS11B protein. (region surrounding Asp88)
Formulation:	Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	1.0 mg/ml
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 36 kDa
Gene Name:	ZFP36 ring finger protein-like 1
Database Link:	<a href="#">Entrez Gene 677 Human Q07352</a>




[View online »](#)

**Background:**


TIS11B (tetradecanoyl phorbol acetate-inducible sequence 11b), also known as ZFP36L1, BRF1 (butyrate response factor 1), ERF1 (EGF response factor 1), cMG1, Berg36 or RNF162B, is a member of the tristetraprolin family. Tristetraprolin (TTP), or TIS11, is a zinc-binding protein encoded by the immediate-early response gene, Zfp-36. TIS11B, a relative of TTP, localizes to the nucleus and may function as a transcription factor involved in regulating the growth factor response. It is an evolutionarily conserved protein containing two C3H1-type zinc fingers and a repeating Cys-His motif. TIS11B is an mRNA binding protein and is known to interact with the 3'-untranslated region of VEGF mRNA, thereby decreasing its stability. This suggests that TIS11B is a potential target in antiangiogenic therapy. In addition, TIS11B may also be an important regulator of myogenesis, as its expression is upregulated during murine myoblast differentiation.

**Synonyms:**

Butyrate response factor 1, BRF1, BERG36, RNF162B, ERF1, ERF-1

**Product images:**


Western blot (WB) analysis of TIS11B antibody (Cat.-No.: [AP01366PU-N]) in extracts from Jurkat cells.



Immunohistochemistry (IHC) analysis of TIS11B antibody (Cat.-No.: [AP01366PU-N]) in paraffin-embedded human brain tissue