

Product datasheet for **AP01427PU-M**

Delta Opioid Receptor (OPRD1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western Blot: 1/500-1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	DOR-1 antibody detects endogenous levels of DOR-1 protein. (region surrounding Phe345)
Formulation:	Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	1.0 mg/ml
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 36 kDa
Gene Name:	opioid receptor delta 1
Database Link:	Entrez Gene 4985 Human P41143

[View online »](#)

Background:

The opioid receptors are G-protein coupled, seven-transmembrane domain receptors for enkephalins, endorphins, and dynorphins. Three different opioid receptor subtypes (kappa , delta, and mu) were first identified by their different selectivities for various naturally occurring alkaloid opioid ligands, and subsequently confirmed by molecular cloning. The amino acid sequences of the opioid receptor subtypes are ~70% homologous, and are similar to somatostatin receptors (SSTRs) showing ~40 % homology with SSTR1. G-protein binding is thought to occur at the third intracellular loop of the opioid receptors, which is also the location of consensus sequences for phosphorylation of the receptor. Interestingly, the genes encoding the specific receptor subtypes are found on different chromosomes in both the human and mouse genomes.

Synonyms:


OPRD1, OPRD, DOR-1, D-OR-1

Protein Families:

Druggable Genome, Transmembrane

Protein Pathways:

Neuroactive ligand-receptor interaction

Product images:


Western blot (WB) analysis of DOR-1 antibody (Cat.-No.: [AP01427PU-N]) in extracts from NIH-3T3 cells.