

Product datasheet for **AP01472PU-N**


FOXB2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Western Blot: 1/500-1/1000.
Reactivity:	Human, Mouse
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	FoxB2 antibody detects endogenous levels of FoxB2 protein. (region surrounding Gln285)
Formulation:	Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	1.0 mg/ml
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 46 kDa
Gene Name:	forkhead box B2
Database Link:	Entrez Gene 14240 Mouse Entrez Gene 442425 Human Q5VYV0
Background:	FOXB2 may act as a transcription factor.
Synonyms:	Forkhead box protein B2, RP11-159H20.4



Product images:



Western blot (WB) analysis of FoxB2 antibody (Cat.-No.: AP01472PU-N) in extracts from HeLa cell