

Product datasheet for **AP01492PU-S**

Oct4 (POU5F1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	Western Blot: 1/500-1/1000.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	OCT3 antibody detects endogenous levels of OCT3 protein. (region surrounding Gln306)
Formulation:	Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	1.0 mg/ml
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 61 kDa
Gene Name:	POU class 5 homeobox 1
Database Link:	Entrez Gene 5460 Human Q01860




[View online »](#)

Background:

Organic cation transporters (OCT) are expressed in the plasma membrane of epithelial cells from a wide range of tissues, where they function in the elimination of endogenous amines, cationic drugs and other xenobiotics. The structure of OCTs consists of a 12-transmembrane-domain structure and a large extracellular hydrophilic loop. In humans, OCT1 is primarily expressed in the liver while OCT2 is expressed in the kidney. OCT3 is expressed in the placenta, skeletal muscle, prostate, aorta and liver. OCT3, also known as extraneuronal monoamine transporter, is widely expressed in different regions of the brain including the hippocampus, cerebellum and cerebral cortex. OCT3 mediates the uptake of several neuroactive agents, including dopamine, and may play an important role in the disposition of neurotransmitters and cationic neurotoxins in the brain. The genes encoding human OCT1–3 map to a conserved cluster at chromosome 6q26-q27.

Synonyms:

POU5F1, OCT3, OCT4, OTF3, POU domain, class 5, transcription factor 1, Octamer-binding transcription factor 3, Oct-3, Oct-4

Product images:

Western blot (WB) analysis of OCT3 antibody (Cat.-No.: [AP01492PU-N]) in extracts from HUVEC cells.