

Product datasheet for **AP01506PU-S**

ADCY4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Western Blot: 1/500-1/1000.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Specificity:	A Cyclase IV antibody detects endogenous levels of A Cyclase IV protein. (region surrounding Leu226)
Formulation:	Phosphate buffered saline (PBS) with 0.05% sodium azide, approx. pH 7.2. State: Aff - Purified State: Liquid purified Ig fraction (> 95% pure by SDS-PAGE)
Concentration:	1.0 mg/ml
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	~ 120 kDa
Gene Name:	adenylate cyclase 4
Database Link:	Entrez Gene 196883 Human Q8NFM4




[View online »](#)

Background:

Adenylyl cyclases function to convert ATP to cyclic AMP in response to activation by a variety of hormones, neurotransmitters and other regulatory molecules. Cyclic AMP, in turn, activates several other target molecules to control a broad range of diverse phenomena such as metabolism, gene transcription and memory. Adenylyl cyclases respond to receptor-initiated signals, mediated by the Gs and Gi heterotrimeric G proteins. The binding of an agonist to a Gs-coupled receptor catalyzes the exchange of GDP (bound to G α s) for GTP, the dissociation of GTP-G α s from G β and G α s-mediated activation of adenylyl cyclase. Adenylyl cyclase IV (AC IV) and IX mRNA are expressed in all kidney nephron segments. AC IV exhibits moderate staining in type II and type IV fibrocytes in rat cochlea and immunoreactivity is also observed in type I fibrocytes. Activation of the D2 dopaminergic and m4 muscarine receptors inhibits the activity of adenylyl cyclase isozymes I, V, VI and VIII, whereas type II, IV and VII are stimulated and type III is not affected.

Synonyms:

Adenylate cyclase type IV, Adenylyl cyclase 4, ADCY4

Product images:

Western blot (WB) analysis of A Cyclase IV antibody (Cat.-No.: [AP01506PU-N]) in extracts from Jurkat cells.