

Product datasheet for **AP02005SU-N**

Asymmetric Dimethylarginine (ADMA) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, R, WB
Recommended Dilution:	RIA. ELISA. Western blot. Immunocytochemistry (1/100). Immunohistochemistry on Cryo Sections. Immunohistochemistry on Paraffin Sections: 1/100.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Asymmetric Dimethylarginine (ADMA) KLH-conjugated
Specificity:	Using conjugated ADMA-KLH, antibody specificity was performed with an ELISA test by competition experiments with the following compounds: Compound : Cross-reactivity ratio(α) ADMA-BSA: 1 SDMA-BSA: 1 / > 100 NDMA-BSA: 1 / > 100 L-Arginine-BSA: 1 / > 5000 (α) ADMA-BSA concentration/unconjugated or conjugated amino acids concentration at half displacement
Formulation:	State: Serum State: Lyophilized Serum.
Reconstitution Method:	Restore in 100 μ l aqua bidest.
Conjugation:	Unconjugated



Storage:	Store lyophilized at 2–8°C for 6 months or at –20°C long term. After reconstitution store the antibody undiluted at 2–8°C for one month or (in aliquots) at –20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Asymetric Dimethyl Arginine
Note:	CAS 65005–57–4