

Product datasheet for **AP02011SU-N**

GM CSF (CSF2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, R
Recommended Dilution:	RIA. Immunohistochemistry.
Reactivity:	Bovine
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Recombinant bovine Granulocyte Macrophage Colony-stimulating Factor
Specificity:	This antibody detects Granulocyte Macrophage Colony-stimulating Factor.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	colony stimulating factor 2
Database Link:	Entrez Gene 1437 Human P04141
Synonyms:	CSF2, GMCSF, Sargramostim, Molgramostin
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Fc epsilon RI signaling pathway, Hematopoietic cell lineage, Jak-STAT signaling pathway, Natural killer cell mediated cytotoxicity, T cell receptor signaling pathway



[View online »](#)