

Product datasheet for **AP02013SU-S**

Growth Hormone (GH1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, R, WB
Recommended Dilution:	RIA (1:60000). Immunohistochemistry (1:1000). Western blot (1:1000).
Reactivity:	Bovine
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Porcine Growth Hormon
Specificity:	This antibody detects Growth Hormon. Shows no cross reactivity with other anterior pituitary hormones.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	growth hormone 1
Database Link:	P01241
Synonyms:	Pituitary growth hormone, Growth hormone 1, HGH



[View online »](#)