

## Product datasheet for **AP02019SU-S**

### **TMSB10 (1-4) Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:500).
Reactivity:	Bovine
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic bovine Thymosin b9 (aa 1-4) KLH conjugated
Specificity:	This antibody detects synthetic Thymosin b9 (aa 1-4) and Thymosin b9.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Database Link:	<a href="#">Entrez Gene 282385 Bovine P21752</a>
Synonyms:	MIG12; PTMB10; TB10; THYB10
Note:	LocusID 282385



[View online »](#)