

Product datasheet for **AP02020SU-S**

TMSB10 (33-41) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:500).
Reactivity:	Bovine
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic bovine Thymosin b9 (aa 33-41) KLH conjugated
Specificity:	This antibody detects Thymosin beta9.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Database Link:	Entrez Gene 282385 Bovine P21752
Synonyms:	MIG12; PTMB10; TB10; THYB10
Note:	LocusID 282385



[View online »](#)