


Product datasheet for **AP02023SU-S**

Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Recommended Dilution:	ELISA, WB
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	N ^ε -(Carboxymethyl)lysine KLH-conjugated
Purification:	Lyophilized serum Serum
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



Product images:

Western Blot analysis of CML-modification in protein extracts of Mouse and Rat tissues. 20 g total protein was separated by SDS-PAGE and immunoblotted with AP02023SU (1/20000) at 4C overnight, followed by a Goat anti Rabbit antibody coupled to HRP (1/10000)