

Product datasheet for **AP02026SU-N**

Eosinophil derived neurotoxin (RNASE2) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1/1000).
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Human Eosinophil-derived Neurotoxin KLH-conjugated
Specificity:	This antibody detects Eosinophil-derived Neurotoxin.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	ribonuclease A family member 2
Database Link:	P10153
Synonyms:	EDN, RNS2, Ribonuclease 2



[View online »](#)