

Product datasheet for **AP02045SU-N**

ANP (NPPA) (ANP 1-28) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human ANP (aa 1-28) poly Lysin conjugated
Specificity:	This antibody detects Atrial Natriuretic Peptide/ANP (aa 1-28).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	natriuretic peptide A
Database Link:	Entrez Gene 4878 Human P01160
Background:	Atrial natriuretic factor (ANP) is a potent vasoactive substance which is thought to play a key role in cardiovascular homeostasis and has a cGMP stimulating activity. There are two different prepronatriodilatin alleles. One codes for 2 Arginine residues at the C terminus that are cleaved to form the mature peptide, while the other ends in a termination codon immediately after the last codon of the mature peptide. ANP belongs to the natriuretic peptide family.
Synonyms:	Atrial Natriuretic Peptide, Prepronatriodilatin, NPPA, ANP, PND
Note:	LocusID 4878



[View online »](#)