

## Product datasheet for **AP02050SU-N**

### Bone Sialoprotein (IBSP) (108-122) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA (1:1000). Western blot (1:250). Immunocytochemistry (1:100).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Bone Sialoprotein (aa 108-122)
Specificity:	This antibody detects Bone Sialoprotein (aa 108-122).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqy bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	integrin binding sialoprotein
Database Link:	<a href="#">Entrez Gene 3381 Human P21815</a>
Background:	The protein encoded by the Bone Sialoprotein gene is a major structural protein of the bone matrix. It is synthesized by skeletal-associated cell types, including hypertrophic chondrocytes, osteoblasts, osteocytes, and osteoclasts, and it constitutes approximately 12% of the noncollagenous proteins in human bone. The only extraskelatal site of its synthesis is the trophoblast. This protein binds to calcium and hydroxyapatite via its acidic amino acid clusters, and mediates cell attachment through an RGD sequence that recognizes the vitronectin receptor.
Synonyms:	Bone sialoprotein II, BSP II, IBSP, BNSP
Note:	LocusID 3381



[View online »](#)