

## Product datasheet for **AP02065SU-N**

### Calcitonin (CALCA) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	RIA (1/15,000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Calcitonin, BSA-conjugated
Specificity:	This antibody detects Calcitonin. There were no cross reactivities obtained with Human Katalcalcin and Human Calcitonin Gene related Peptide.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	calcitonin related polypeptide alpha
Database Link:	<a href="#">Entrez Gene 796 Human P06881</a>
Background:	Calcitonin is a 32 amino acid peptide hormone synthesized by the parafollicular cells of the thyroid. It causes a rapid, but short lived, reduction in serum calcium and phosphate by promoting the incorporation of those ions in the bones. This effect is opposite to that of parathyroid hormone. Staining for calcitonin may be used for the identification of a spectrum of C cell proliferative abnormalities ranging from C cell hyperplasia to invasive tumors. Staining for calcitonin in medullary carcinoma of the thyroid produces a fine granular pattern in the cytoplasm. Amyloid deposits within the tumor may also exhibit varying degrees of calcitonin activity.
Synonyms:	CALCA, CALC1



[View online »](#)

Note: LocusID 796