

Product datasheet for **AP02076SU-S**

Macrophage inflammatory protein 5 (CCL15) (18-31) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R, WB
Recommended Dilution:	ELISA (1:500). RIA (1:1000). Western blot (1:500).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Chemokine CC-2 (aa 18-31) poly Lysinconjugated
Specificity:	This antibody detects Chemokine CC-2 (aa 18-31)
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	C-C motif chemokine ligand 15
Database Link:	Entrez Gene 6359 Human Q16663
Background:	Macrophage inflammatory protein 5 (CCL15) is a member of the human beta-chemokine family and is chemotactic for T cells and monocytes but not neutrophils, eosinophils, or B cells. It induces changes in intracellular calcium concentration in monocytes and is thought to act through the CCR1 receptor but also binds to CCR3.
Synonyms:	C-C motif chemokine 15, Small-inducible cytokine A15, Chemokine CC-2, NCC-3, NCC3, MIP-5, MIP-1 delta, Leukotactin-1, Mrp-2b, SCYA15
Note:	LocusID 6359



[View online »](#)