

Product datasheet for AP02092SU-S

Endothelin 1 (EDN1) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA, IHC, R, WB

Recommended Dilution: RIA.

EIA.

Western blot.

Immunohistochemistry on Paraffin Sections.

Reactivity: Bovine, Human, Porcine

Host: Rabbit

Clonality: Polyclonal

Immunogen: Human Endothelin-1

Specificity: This antibody detects Endothelin-1.

Formulation: State: Serum

State: Lyophilized Serum without additives

Reconstitution Method: Restore in aqua bidest to initial volume.

Conjugation: Unconjugated

Storage: Prior to reconstitution store at 2-8°C.

Following reconstitution store undiluted at 2-8°C for one month

or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: endothelin 1

Database Link: Entrez Gene 1906 Human

P05305

Background: Endothelins (ET) show potent constrictor activity in vascular and non-vascular smooth

muscle. This family of 21-amino acid peptides exists in at least three isoforms - ET-1, ET-2, and ET-3, and is produced in endothelial and epithelial cells. ET's can mediate biological effects in cells and tissues, and have been shown to bind to an ET receptor in the lung, kidney, heart, and liver. Endothelin 1 is expressed in lung, placental stem villi vessels and in

cultured placental vascular smooth muscle cells.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com





Synonyms: PPET1, EDN1, ET1, HDLCQ7, EDN1