

Product datasheet for **AP02132SU-N**

Osteocalcin (BGLAP) (1-22) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R
Recommended Dilution:	RIA (1/2000). ELISA (1/5000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Osteocalcin (aa 1-22)
Specificity:	This antibody detects Osteocalcin (aa 1-22). There were no cross reactivities obtained with Human Osteopontin, Human Osteonectin, and Human Bone Sialoprotein.
Formulation:	State: Serum State: Lyophilised Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	bone gamma-carboxyglutamate protein
Database Link:	Entrez Gene 632 Human P02818
Synonyms:	Bone Gla protein, BGLAP
Note:	LocusID 632



[View online »](#)