

## Product datasheet for **AP02132SU-S**

### Osteocalcin (BGLAP) (1-22) Rabbit Polyclonal Antibody

#### Product data:

|                        |   |
|------------------------|---|
| Product Type:          | Primary Antibodies  |
| Applications:          | ELISA, R  |
| Recommended Dilution:  | <b>RIA</b> (1/2000).<br><b>ELISA</b> (1/5000).  |
| Reactivity:            | Human   |
| Host:                  | Rabbit  |
| Clonality:             | Polyclonal  |
| Immunogen:             | Synthetic Human Osteocalcin (aa 1-22)   |
| Specificity:           | This antibody detects Osteocalcin (aa 1-22).<br>There were no cross reactivities obtained with Human Osteopontin, Human Osteonectin, and Human Bone Sialoprotein. |
| Formulation:           | State: Serum<br>State: Lyophilised Serum  |
| Reconstitution Method: | Restore in aqua bidest to initial volume.   |
| Conjugation:           | Unconjugated  |
| Storage:               | Store lyophilized at 2-8°C and reconstituted at -20°C.<br>Avoid repeated freezing and thawing.  |
| Stability:             | Shelf life: One year from despatch.   |
| Gene Name:             | bone gamma-carboxyglutamate protein   |
| Database Link:         | <a href="#">Entrez Gene 632 Human P02818</a>  |
| Synonyms:              | Bone Gla protein, BGLAP   |
| Note:                  | LocusID 632   |



[View online »](#)