

Product datasheet for **AP02134SU-N**

SPARC (272-286) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1/3000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Osteonectin (aa 272-286)
Specificity:	This antibody detects Synthetic Osteonectin (aa 272-286), and recombinant Osteonectin. There were no cross reactivities obtained with human Osteopontin, and human Bone Sialoprotein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	secreted protein acidic and cysteine rich
Database Link:	Entrez Gene 6678 Human P09486
Synonyms:	ON, Basement-membrane protein 40, BM-40
Note:	LocusID 6678



[View online »](#)