

## Product datasheet for **AP02143PU-S**

### Parathyroid Hormone (PTH) (1-34) Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	<b>RIA</b> (1/30000).
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Synthetic Human PTH (aa 1-34), bTG- conjugated
Specificity:	This antibody detects PTH (aa 12-17; 1-34; 1-38; 1-84).
Formulation:	PBS State: Ig Fraction State: Lyophilized purified IgG fraction
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Purification:	Protein G Chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	parathyroid hormone
Database Link:	<a href="#">Entrez Gene 5741 Human P01270</a>
Synonyms:	Parathormone, Parathyrin
Note:	LocusID 5741



[View online »](#)