

## Product datasheet for **AP02143SU-S**

### Parathyroid Hormone (PTH) (1-34) Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R
Recommended Dilution:	<b>ILMA.</b> <b>ELISA:</b> 1/1000 (See Ref.1)
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Synthetic Human PTH (aa 1-34), bTG- conjugated
Specificity:	This antibody detects PTH (aa 12-17, 1-34, 1-38, 1-84).
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8 °C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	parathyroid hormone
Database Link:	<a href="#">Entrez Gene 5741 Human P01270</a>
Synonyms:	Parathormone, Parathyrin
Note:	LocusID 5741



[View online »](#)