

Product datasheet for **AP02151SU-S**

Parathyroid Hormone (PTH) (28-48) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	RIA (1/30,000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human PTH (aa 28-48)
Specificity:	This antibody detects PTH (aa 28-48).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	parathyroid hormone
Database Link:	Entrez Gene 5741 Human P01270
Synonyms:	Parathormone, Parathyrin



[View online »](#)