

## Product datasheet for **AP02153SU-S**

### Parathyroid Hormone (PTH) (53-84) Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R
Recommended Dilution:	RIA (1:10,000). ELISA (1:5,000).
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Synthetic human PTH (aa 53-84) bTG-conjugated
Specificity:	This antibody detects PTH (aa 53-68; 53-84; 1-84).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	parathyroid hormone
Database Link:	<a href="#">Entrez Gene 5741 Human P01270</a>
Synonyms:	Parathormone, Parathyrin
Note:	LocusID 5741



[View online »](#)