

Product datasheet for **AP02154SU-N**

Parathyroid Hormone (PTH) (53-84) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	RIA (1:5.000).
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human PTH (aa 53-84) bTG-conjugated
Specificity:	This antibody detects PTH (aa 53-68; 53-84; 1-84).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	parathyroid hormone
Database Link:	Entrez Gene 5741 Human P01270
Synonyms:	Parathormone, Parathyrin
Note:	LocusID 5741



[View online »](#)