

## Product datasheet for **AP02160PU-S**

### Parathyroid hormone related protein (PTHrP) (1-34) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	<b>RIA.</b>
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic Human PTHrP (aa 1-34) KLH conjugated
Specificity:	This antibody detects Human PTHrP (aa 1-34). Other species not tested.
Formulation:	PBS, pH 7.4 State: Purified State: Lyophilized purified IgG fraction from Serum
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml
Purification:	Affinity Chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	parathyroid hormone-like hormone
Database Link:	<a href="#">Entrez Gene 5744 Human P12272</a>
Synonyms:	PTH-rP, Parathyroid hormone related protein, PTHLH, PLP, Parathyroid hormone-like protein



[View online »](#)