

## Product datasheet for **AP02161PU-N**

### Parathyroid hormone related protein (PTHrP) (53-86) Goat Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	<b>RIA</b> (Ref.1)
Reactivity:	Human
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic Human PTHrP (aa 53-86) KLH-conjugated
Specificity:	This antibody detects PTHrP (aa 53-86; 1-86).
Formulation:	PBS, pH 7.4 State: Purified State: Lyophilized purified IgG fraction from Serum
Reconstitution Method:	Restore in 100 µl aqua bidest
Purification:	Protein G Chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	parathyroid hormone-like hormone
Database Link:	<a href="#">Entrez Gene 5744 Human P12272</a>
Synonyms:	PTH-rP, Parathyroid hormone related protein, PTHLH, PLP, Parathyroid hormone-like protein
Note:	LocusID 5744



[View online »](#)