

Product datasheet for **AP02161SU-N**

Parathyroid hormone related protein (PTHrP) (53-86) Goat Polyclonal Antibody

Product data:

| | |
|------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | R |
| Recommended Dilution: | RIA (1:2000). |
| Reactivity: | Human |
| Host: | Goat |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic human PTHrP (aa 53-86) KLH-conjugated |
| Specificity: | This antibody detects PTHrP (aa 53-86; 1-86). |
| Formulation: | State: Serum State: Lyophilized serum |
| Reconstitution Method: | Restore in aqua bidest to initial volume. |
| Conjugation: | Unconjugated |
| Storage: | Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: One year from despatch. |
| Gene Name: | parathyroid hormone-like hormone |
| Database Link: | Entrez Gene 5744 Human P12272 |
| Synonyms: | PTH-rP, Parathyroid hormone related protein, PTHLH, PLP, Parathyroid hormone-like protein |
| Note: | Reference no.: LocusID 5744 |



[View online »](#)