

Product datasheet for **AP02161SU-S**

Parathyroid hormone related protein (PTHrP) (53-86) Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	RIA (1:2000).
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Synthetic human PTHrP (aa 53-86) KLH-conjugated
Specificity:	This antibody detects PTHrP (aa 53-86; 1-86).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	parathyroid hormone-like hormone
Database Link:	Entrez Gene 5744 Human P12272
Synonyms:	PTH-rP, Parathyroid hormone related protein, PTHLH, PLP, Parathyroid hormone-like protein
Note:	Reference no.: LocusID 5744



[View online »](#)