

Product datasheet for **AP02163SU-N**

SERPINB2 Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA (1:5000). Western blot (1:100). Immunohistochemistry on paraffin (5 - 10 µg/ml) and frozen sections.
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Purified human PAI-2
Specificity:	This antibody detects PAI 2. There were no cross reactivities obtained with human tPA ; uPA; and PAI-1.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	serpin family B member 2
Database Link:	Entrez Gene 5055 Human P05120
Synonyms:	Serpin B2, PAI-2, Plasminogen activator inhibitor 2, PLANH2
Note:	Reference no.: LocusID 5055



[View online »](#)