

## Product datasheet for **AP02184SU-N**

### Guanylin (GUCA2A) (4-16) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, R, WB
Recommended Dilution:	ELISA. RIA. Western blot.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human pro-Guanylin (aa 4-16)
Specificity:	This antibody detects synthetic pro-Guanylin (aa 4-16).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	guanylate cyclase activator 2A
Database Link:	<a href="#">Entrez Gene 2980 Human Q02747</a>
Synonyms:	Guanylate cyclase activator 2A, GUCA2A, GUCA2
Note:	Reference no.: LocusID 2980



[View online »](#)