

## Product datasheet for **AP02186SU-N**

### PRAP1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/10.000. Western blot: 1/15.000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic human Proline-rich acidic protein 1 (PRAP1) Peptide (aa 58-69), polyLysine-conjugated
Specificity:	This antibody detects Proline-rich acidic protein 1 (PRAP1) Peptide (aa 58-69).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	proline-rich acidic protein 1
Database Link:	<a href="#">Entrez Gene 118471 Human Q96NZ9</a>
Background:	PRAP1 (Proline rich acidic protein 1) is specifically expressed in the intestine and uterus. It may play a role in pregnant uterus, and in digestion and absorption of nutrients in the gastrointestinal tract.
Synonyms:	MGC126792; PRO1195; RP11-122K13.6; UPA



[View online »](#)