

Product datasheet for **AP02188SU-N**

NPPC (proCNP 1-19) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1:1.000).
Reactivity:	Human
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Synthetic human pro-CNP (aa 1-19)
Specificity:	This antibody detects pro CNP (aa 1-19).
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	natriuretic peptide C
Database Link:	Entrez Gene 4880 Human P23582
Synonyms:	NPPC, CNP2
Note:	Reference no.: LocusID 4880



[View online »](#)