

Product datasheet for **AP02189SU-S**

NPPC (proCNP 30-50) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1/1000).
Reactivity:	Human
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Synthetic Human pro-CNP (aa 30-50)
Specificity:	This antibody detects synthetic Human pro-CNP (aa 30-50).
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	natriuretic peptide C
Database Link:	Entrez Gene 4880 Human P23582
Synonyms:	NPPC, CNP2



[View online »](#)