

## Product datasheet for **AP02190SU-S**

### **NPPC (proCNP 51-80) Sheep Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (1/1,000).
Reactivity:	Human
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Synthetic Human pro-CNP (aa 51-79)
Specificity:	This antibody detects synthetic pro-CNP (aa 51-79).
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	natriuretic peptide C
Database Link:	<a href="#">Entrez Gene 4880 Human P23582</a>
Synonyms:	NPPC, CNP2
Note:	Reference no.: LocusID 4880



[View online »](#)