

Product datasheet for **AP02200SU-N**

Relaxin 2 (RLN2) (175-185) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Relaxin (aa 175-185) KLH conjugated
Specificity:	This antibody detects Relaxin (aa175-185) (= Relaxin A chain, aa14-24).
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	relaxin 2
Database Link:	Entrez Gene 6019 Human P04090
Synonyms:	Relaxin A chain, Relaxin B chain, RLN1, RLN2, RLN, Prorelaxin H1, Prorelaxin H2



[View online »](#)