

## Product datasheet for **AP02201SU-S**

### Relaxin 2 (RLN2) (28-34) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic Human Relaxin (aa 28-34) KLH conjugated
Specificity:	This antibody detects Relaxin (aa28-34) (= Relaxin B chain, aa4-10).
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in aqua bidest to initial volume.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C and reconstituted at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	relaxin 2
Database Link:	<a href="#">Entrez Gene 6019 Human P04090</a>
Synonyms:	Relaxin A chain, Relaxin B chain, RLN1, RLN2, RLN, Prorelaxin H1, Prorelaxin H2



[View online »](#)