

Product datasheet for **AP02207PU-S**

Resistin (RETN) (30-49) Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA (3 µg/ml).
Reactivity:	Human
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Synthetic human Resistin (aa 30-49) bTG- conjugated
Specificity:	This antibody detects Resistin.
Formulation:	50 mM TRIS pH 7,4 State: Purified State: Lyophilized antibody
Reconstitution Method:	Restore in aqua bidest to 1 mg/ml.
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Gene Name:	resistin
Database Link:	Entrez Gene 56729 Human Q9HD89
Synonyms:	RETN, FIZZ3, HXCP1, RSTN, ADSF



[View online »](#)