

Product datasheet for **AP02209SU-N**

Resistin (RETN) (33-90) Rabbit Polyclonal Antibody

Product data:

| | |
|------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | ELISA |
| Recommended Dilution: | ELISA. |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic human Resistin (aa 33-90) bTG-conjugated |
| Specificity: | This antibody detects Resistin. |
| Formulation: | State: Serum State: Lyophilized serum |
| Reconstitution Method: | Restore in aqua bidest to initial volume. |
| Conjugation: | Unconjugated |
| Storage: | Store lyophilized at 2 - 8 °C and reconstituted at -20 °C. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: One year from despatch. |
| Gene Name: | resistin |
| Database Link: | Entrez Gene 56729 Human Q9HD89 |
| Synonyms: | RETN, FIZZ3, HXCP1, RSTN, ADSF |



[View online »](#)